

# UNIT COUNTING GUIDE




**CLASPS**  
 Akers, 1 unit each  
 RPI, RPA, RPR, and Equipoise  
 1½ Units each



**HORIZONTAL FINISH LINES**  
 3 to 6 teeth including retention  
 1½ units, 7 or more teeth or acrylic  
 palate, 2 units.




**CLASPS WITH LUGS**  
 1½ Units each



**KENNEDY BAR**  
 ¼ unit per tooth passed



**RESTS**  
 Plain rest ¼ unit  
 Individual res with strut ½ unit



**APRON OR LINGUAL PLATE**  
 3 units stippled  
 4 units smooth



**LINGUAL BAR**  
 2 units



**HORSESHOE**  
 For each tooth passed, ¼ unit  
 2 units maximum



**STRENGTHENING BARS**  
 4 units



**CIRCULAR BAR**  
 Double palatal bars 2 units  
 for each bar



**DISTAL RETENTIONS**  
 With vertical finish line 1 unit



**PALATAL BARS**  
 3 units



**HORIZONTAL FINISH LINE**  
 Up to 2 teeth including retention  
 1 unit



**LOWER CAST BASE**  
 4 units



**UPPER CAST BASE**  
 5½ units. With mesh post dam  
 6½ units

## Turn Around Times

Frames = Next Day  
 Frames with Acetal Resin Clasp = Next Day + 1  
 Frames with Metal Teeth, Tube Teeth, Onlays, or Bite Openers = 4 Days or less  
 Flexibles = 4 Days